



SophistiCart

Stop bacteria!

The zero downtime cart is also anti-bacterial coated



Automatically adjusted



Best user experience



Modular



Zero downtime



Lightweight



Anti-bacterial coating




Sophisticart is an innovative-tech cart. Automatically adjusted, modular, ergonomic and lightweight. It is the first anti-bacterial coated cart in the world and the best user experience.

No more downtime. No more bulky carts - welcome the next generation of carts.



See the game changer advantages:



EASY TO INSTALL MONITOR

VESA tilt (up/down left/right) that slides into monitor's arm bay, and allows easy installation of any AIO pc.



LARGE WORKSPACE

50cm x 50cm, anti bacterial ABS material for safe MEDprep. Easy to clean and implement infection control.



SLIDE OUT KEYBOARD DRAWER

With left/right mouse pad and UV LIGHT (charging indication) to keep keyboard and mouse germ free.



FLEXIBLE

Convert any basic cart into a medication dispensing cart by adding a 4 or 5 tier medication cabinet with electric pin code lock. Accessories are easily attachable to the cart



EASY HEIGHT ADJUSTMENT BY ELECTRIC LIFT

Suitable for a load capacity of 100 kg



GREEN BATTERY

For 13, and up to 24 working hours. Fast charging, long life and lightweight. Easy maintenance, simply replace the compartment in minutes in the event of failure (no downtime). Clear LCD battery status indication.

See the game changer details:



TECHNICAL SPECIFICATION

- Net Weight: 25 Kg
- Base Size: 60cm x 50cm
- Working surface Height Range: 81cm-126cm
- Height adjustment Range: 45cm
- Work Surface: 50cm x 50cm , anti bacterial
- Monitor arm with VESA tilt- up to 24" monitor
- Casters: 5" (12.7 cm), 2 locking
- Lift Type: Electronic. 1000N (about 100kg)
- Display Max Weight: 16Kg
- Battery: Lifepo4 300W or 600w, 3000 life cycles
- Power input: 110-220v
- Power output: 12v/19v
- Germicidal UV light source within keyboard compartment
- Warranty: 2 years

ORDERING INFORMATION

- Easy - 300 Electric lift, Life powered cart 300w
- Easy - 600 Electric lift, Life powered cart 600w
- Easy - 300M Medication cart, Electric lift, Life powered cart 300w
- Easy - 600M Medication cart, Electric lift, Life powered cart 600w



www.mennenmedical.com

Israel: Mennen Medical Ltd
6 Ha-Klshon Street, Yavne
P.O. Box 102, Rehovot 7610002, Israel
Phone: +972-8-9323333
Fax: +972-8-9328510

United States: Mennen Medical Corp
290 Andrews Road
Feasterville-Treose, PA 19053-3480, USA
Phone: +1-215-259-1020
Fax: +1-215-357-2010

United Kingdom: Charter Kontron Ltd
Unit 18 Avant Business Centre
21 Denbigh Road Milton Keynes MK1 1DT, UK
Phone: +44-1908-646070
Fax: +44-1908-646030